

KOMET® Sonic Tips

The KOMET® SonicLine sonic-handpiece system now includes special sonic tips for the opening of fissures, for the finishing of veneer preparations, and for crown preparation.

Operated at low oscillation amplitude for minimally invasive fissure-opening procedures, the KOMET® SF849 sonic tip features a small-diameter tip that allows the operator to detect hidden caries, to remove fissure caries, and to prepare fissures for sealant applications.

The KOMET® SF8850 sonic tip for veneer-preparation procedures is operated at low oscillation amplitude, yielding the exceptionally fine surface necessary to provide a tight margin for the veneer prep.

The KOMET® SonicLine crown-preparation tips are offered in six configurations to produce precisely positioned and finished prosthetic margins and to enhance interproximal finishing. The working ends of two additional configurations, specially indicated for interproximal trimming, are half-coated lengthwise for applications on mesial and distal surfaces.

All of KOMET's sonic tips are designed exclusively for use with the KOMET® sonic handpiece (SF1LM).

SonicLine Tips

Fissures



SF 849.000.009 \$125.00

Veneers



SF 8850.000.016 \$125.00

Crown Preparation and Finishing ●



SF 979.000.014 \$125.00
SF 8979.000.014 ● \$125.00



SF 847KR.000.016 \$125.00



SF 979.000.016 \$125.00
SF 8979.000.016 ● \$125.00



SF 8878KD.000.018 ● \$125.00



SF 862.000.014 \$125.00



SF 8878KM.000.018 ● \$125.00

● = Fine Grit

"KOMET's SonicLine sonic tips have been a great addition to my armamentarium. They allow me to perform dentistry safely, precisely, and with unmatched reward. The sonic tips simplify the preparation of subgingival finish lines for crowns and veneers, and they offer the added benefit of same-day, blood-free, truly accurate impressions."

Alessandro Geminiani DDS, Rochester, NY



Fissures



Veneers



Crown Preparation and Finishing

