

S-Diamonds®

Quick and effective
preparation



- Increased efficiency
- Combination of a structured blank and coarse grit
- Quick and effective substance removal saves time
- Measurably higher material reduction when using structured instruments
- Patented design

A world-wide leader
in rotary instruments

888-566-3887
www.komet-usa.com




Feel the impact


Product characteristics and advantages

The primary preparation of crown cores is one of the most time consuming steps during treatment. The use of S-Diamonds® is highly recommended for increased efficiency and gentle treatment. The combination of a structured blank and coarse grit allows quick and effective substance removal and improved cooling. The faceted structure of the blank reduces clogging and the generation of heat. All in all, the operation becomes notably more effective.

The time saved during treatment is equally beneficial to you and your patients!



Material reduction in %



A comparison between conventional diamond instruments and S-Diamond® instruments clearly shows a measurably higher material reduction when using structured instruments.

<ul style="list-style-type: none"> S 6856.FG.014 S 6856.FG.016 S 6856.FG.018 S 6856.FG.021 S 5856.FG.014 S 5856.FG.016 S 5856.FG.018 S 5856.FG.021 	<ul style="list-style-type: none"> S 6882.FG.012 S 6882.FG.014 S 6882.FG.016 S 6847KR.FG.014 S 6847KR.FG.016 S 6847KR.FG.018 S 5847KR.FG.014 S 5847KR.FG.016 S 5847KR.FG.018 	<ul style="list-style-type: none"> S 6836KR.FG.012 S 6836KR.FG.014 S 6836KR.FG.016 S 6837KR.FG.012 S 6837KR.FG.014 S 6837KR.FG.016 S 6878K.FG.012 S 6878K.FG.014 S 6878K.FG.016 S 6878K.FG.018 S 6878K.FG.016 S 6878K.FG.018 S 6878K.FG.016 S 5878K.FG.014 S 5878K.FG.016 S 5878K.FG.018 	<ul style="list-style-type: none"> S 6877.FG.012 S 6878.FG.012 S 6878.FG.014 S 6878.FG.016 S 6879.FG.012 S 6879.FG.014 S 6879.FG.016 S 6879.FG.012 S 6879.FG.014 S 6879.FG.016 S 6886K.FG.018 S 6379.FG.018 S 6379.FG.023 S 5379.FG.018 S 5379.FG.023 	<ul style="list-style-type: none"> S 6801.FG.012 S 6801.FG.014 S 6801.FG.016 S 6801.FG.018 S 6801.FG.023 S 6830L.FG.012 S 6830L.FG.014 S 6830L.FG.016 S 6830L.FG.018 S 6830RL.FG.012 S 6830RL.FG.014 S 6830RL.FG.016 S 6830RL.FG.018 	<ul style="list-style-type: none"> S 6845KR.FG.018 S 6845KR.FG.025 S 6835KR.FG.012 S 6835KR.FG.014
--	---	--	--	---	--

Diamond grit sizes

- coarse 151µm
- super-coarse 181µm

US patent 6,368,107
© KOMET USA LLC 11/08 403506V2

US patent 6,368,107

888-566-3887
www.komet-usa.com



Feel the impact